

Download «Microlife Connected Health+» App for free from Google Play™ (Android) or App Store (iOS) and install it on your smartphone.

How the Bluetooth® Function works

The Bluetooth® function on your device will automatically turn on and will be ready to connect with the «Microlife Connected Health+» App after the device is switched on. Your device will automatically upload the data once it is connected to the smartphone.

10. Malfunctions and actions to take

Description	Symptom/ Possible causes	Solutions
SpO ₂ or pulse rate do not display normally.	<ol style="list-style-type: none"> 1. Finger is not inserted correctly. 2. Patient SpO₂ value is too low to be measured. 3. There is excessive illumination. 	<ol style="list-style-type: none"> 1. Retry inserting the finger. 2. & 3. Measure more times. If you determine the product is working correctly, consult your doctor.
SpO ₂ or pulse rate is shown unstable.	<ol style="list-style-type: none"> 1. Finger might not be inserted deep enough. 2. Excessive patient movement. 	<ol style="list-style-type: none"> 1. Retry inserting the finger. 2. Sit calmly and retry.

Description	Symptom/ Possible causes	Solutions
The device cannot be powered on.	<ol style="list-style-type: none"> 1. No batteries or low battery power. 2. Batteries are not installed correctly. 3. The device may be damaged. 	<ol style="list-style-type: none"> 1. Replace the batteries. 2. Remove and reinstall the batteries. 3. Contact your local Microlife-Customer Service.
The display suddenly switches off.	<ol style="list-style-type: none"> 1. The device is automatically powered off, when no signal was detected after 10±2 seconds. 2. The battery power is too low to operate. 	<ol style="list-style-type: none"> 1. Normal. 2. Replace the batteries.